

FLOATING TOP TABLE

Portable x-ray Generator Floating top table





Suitable for use with our full range of x-ray generators, our floating top table now enables you to accurately and simply get perfect positioning every time.

The radiolucent table top offers movement in every horizontal direction possible. After positioning, release the foot operated lock, and the top locks noiselessly into position.

There is now no need to disturb the patient when the cassette is changed or a grid is used, as both are held in trays underneath the table top.

As the x-ray generator is always automatically centred onto the cassette, beam alignment is assured every time.

- Under surface self centring cassette tray
- Under surface grid tray
- Standard table top 1.6m long
- Optional 2.0m table top
- Suitable for all film and CR cassettes as well as all VXR DR panels

Specification

Code UN1007 Floating top table 1600 x 650mm Floating top table 2000 x 650mm

www.veterinaryxrays.co.uk



Regent Park, Princes Risborough Buckinghamshire. HP27 9LE

- t. 01844 344 911
- f. 01844 344 915
- e. office@veterinaryxrays.co.uk

Products and specifications are subject to change without prior notice. E&OE Copyright Veterinary X-rays 2015