

Floating Top Table

Suitable for use with our full range of x-ray generators, the radiolucent tabletop moves freely in all horizontal directions, allowing precise positioning with minimal patient disturbance.

The DR panel is positioned directly beneath the tabletop and aligned to the generator at all times.

KEY FEATURES



Under surface self centring DR panel tray



Complete with VXR generator support



Silent top braking system



Standard table top 1600x650mm



Optional table top 2000x650mm



Suitable for all VXR DR panels or standard cassettes sizes



FLOATING TOP TABLE SPECIFICATION

Code	UN1007	Floating top table 1600 x 650mm
	UN1008	Floating top table 2000 x 650mm

